



Laboratory Products



**Introducing the
latest addition to our
Platinum Series Columns.**

**Discover how
CHIRALPAK IC is the
solution you've been
looking for.**

New cellulosic chiral selector

Immobilized phase for robust column and full
solvent compatibility

Platinum Series of immobilized chiral columns —
CHIRALPAK® IA™, IB™, and IC™



CHIRAL
TECHNOLOGIES INC

CHIRAL TECHNOLOGIES, INC., West Chester, PA USA


1-800-6CHIRAL TEL 610-594-2325 FAX

www.chiraltech.com




4-Benzoyloxy-2-azetidinone

CHIRALPAK® IC™ 4.6 x 250 mm
Mobile Phase: Hexane:Ethanol:DEA 80:20:0.1




Devrinol

CHIRALPAK® IC™ 4.6 x 250 mm
Mobile Phase: Hexane:Ethyl Acetate:DEA 70:30:0.1




Laudanosine

CHIRALPAK® IC™ 4.6 x 250 mm
Mobile Phase: Methanol:0.1% DEA




Diperodon

CHIRALPAK® IC™ 4.6 x 250 mm
Mobile Phase: Hexane:2-Propanol:DEA 80:20:0.1




Praziquantel

CHIRALPAK® IC™ 4.6 x 250 mm
Mobile Phase: THF



Indoprofen

CHIRALPAK® IC™ 4.6 x 250 mm
Mobile Phase: CH₂Cl₂:MeOH:TFA 98:2:0.1



PART NUMBER	PRODUCT NAME	INTERNAL DIAMETER (mm)	COLUMN LENGTH (mm)	PARTICLE SIZE (µ)	PRODUCT TYPE
83324	CHIRALPAK IC	4.6	150	5	Analytical
83325	CHIRALPAK IC	4.6	250	5	Analytical
83311	CHIRALPAK IC (pack of 3)	4.0	10	5	Guard Cartridges



CHIRAL
TECHNOLOGIES INC

CHIRAL TECHNOLOGIES, INC.
West Chester, PA USA

1-800-6CHIRAL TEL
610-594-2325 FAX

www.chiraltech.com